

Oxygen simplified

EverFlo oxygen concentrator

The smaller, lighter EverFlo oxygen concentrator from Philips Respironics is making a big impression on providers and patients. Its sleek, ergonomic design makes EverFlo easier to store, move, and carry than conventional concentrators. And it is one of the quietest and most energy efficient devices available today. With no filter for patients to maintain, a durable metal oxygen outlet, and fewer moving parts than other concentrators, EverFlo is designed to deliver consistent, trouble-free operation.





Molded handle that's easy to pick up and transport with one or both hands



Metal cannula connector for greater durability



Humidifier bottle holder compatible with all bottle styles



Recessed flow meter reduces the risk of breakage



Internal tubing storage keeps tubing conveniently available for use



Interface panel with on/off switch and alarm indicators

| Product specifications | | | |
|--|--|---|---|
| | Standard version | Low power version | IKK version |
| Power | 230V/50Hz | 230V/50Hz | 230V/50Hz |
| Liter flow | 0.5-5 liters per min | 0.5-5 liters per min | 0.5-5 liters per min |
| Sound level | 43 dBA | 43 dBA | < 40 dBA |
| Oxygen purity | 90% - 96% | 87% - 96% | 87% - 96% |
| Weight | 14 kg (32 lbs) | 14 kg (32 lbs) | 15 kg (33 lbs) |
| Dimensions (HxWxD) | 580mm (22.8"), 380mm (15"), 240mm (9.5") | 580mm (22.8"), 380mm (15"), 240mm (9.5") | 580mm (22.8"), 380mm (15"), 240mm (9.5") |
| OPI (oxygen percentage indicator) alarm levels | Low oxygen: 82%, very low oxygen: < 70% | Low oxygen: 82%, very low oxygen: < 70% | Low oxygen: 82%, very low oxygen: < 70% |
| Required maintenance | Compressor intake filter: 2 yrs. Micro disk filter: as needed | Compressor intake filter: 2 yrs. Micro disk filter: as needed | Compressor intake filter: 2 yrs. Micro disk filter: as needed |
| Power consumption | < 300 W | < 300 W | < 300 W |
| Part numbers | 1020006 Europe 1020009 South Africa 1020011 Italy 1020010 Australia 1020012 Argentina 1020013 South Korea 1020016 Hong Kong 1020017 Switzerland | 1020008 UK, Ireland | 1020007 Germany |

Philips Respironics and EverFlo are trademarks of Koninklijke Philips N.V. and its affiliates. All rights reserved.



©2015 Koninklijke Philips N.V. All rights are reserved.

Philips Healthcare reserves the right to make changes in specifications and/or to discontinue any product at any time without notice or obligation and will not be liable for any consequences resulting from the use of this publication. Philips Respironics Asia Pacific +65 6882 5282 Philips Respironics Australia +61 (2) 9947 0440 Philips Respironics United Kingdom +44 800 1300 845

Philips Respironics 1010 Murry Ridge Lane Murrysville, PA 15668 USA

www.philips.com/respironics

Broudy HT 11/23/15 MCI 4107077 PN 1109098