

# Compressor Nebulizer

Item#

**NEB100**

Comes with Disposable Adult Neb Kit and Filters

**NEB100-R**

Comes with Disposable and Reusable Adult Neb Kits and Filters

Neb kits include nebulizer bottle, t-piece, mouthpiece, 6in corrugated reservoir tube, and 7ft supply tube.

HCPCS# E0570 + A7003 + A7013



Item#

**NEB100-BAG**

Carry bag



Great value

**FAST  
TREATMENT  
TIMES**

## Specifications

**Weight** 2.75 lbs

**Input** 120V/60Hz

**Current** 1.5A

**Medication Capacity** 6ml

**Particle Size** Respirable Fraction 0.3-0.4um 0.5um 80%

**MMAD** 2.0um

**Average Nebulization Rate** 0.3-0.4ml/min

**Noise Level** Below 60dBA

**Compressor Pressure Range** 30-45psi

**Operation Pressure Range** 12-14.5psi (0.85-1.1 bar)


**Liter Flow Range** 8-10LPM

**Operation Temperature Range** 50°F-104°F (10°C-40°C)

**Operation Humidity Range** 10 to 95% RH

**Storage Temperature Range** 14°F-149°F (-10°C-65°C)

**Storage Humidity Range** 10 to 95% RH



Comes in an attractive retail box

